


Product Datasheet

TRIM9 antibody (orb1333195)

Description	Mouse monoclonal antibody to TRIM9
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 284-669 of human TRIM9 (NP_443210) produced in E.coli.
Target	TRIM9
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	61.2 kDa
Uniprot ID	Q9C026
NCBI	NM_052978
Dilution Range	WB 1:1000, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...