


Product Datasheet

XTP4 (MIEN1) antibody (orb1333232)

Description	Mouse monoclonal antibody to XTP4 (MIEN1)
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Full length human recombinant protein of human C17orf37 (NP_115715) produced in E.coli.
Target	MIEN1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	12.2 kDa
Uniprot ID	Q9BRT3
NCBI	NM_032339
Dilution Range	WB 1:2000, IHC 1:150, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...