


Product Datasheet

ZFAND3 antibody (orb1333251)

Description	Mouse monoclonal antibody to ZFAND3
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human ZFAND3(NP_068762) produced in HEK293T cell.
Target	ZFAND3
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	25 kDa
Uniprot ID	Q9H8U3
NCBI	NM_021943
Dilution Range	WB 1:4000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...