


Product Datasheet

SEPTIN1 antibody (orb1333287)

Description	Mouse monoclonal antibody to SEPTIN1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human SEPT1(NP_443070) produced in HEK293T cell.
Target	SEPTIN1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	41.8 kDa
Uniprot ID	Q8WYJ6
NCBI	NM_052838
Dilution Range	WB 1:200~1000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...