

## Product Datasheet

### Munc 13 4 (UNC13D) antibody (orb1333294)

**Description** Mouse monoclonal antibody to Munc 13 4 (UNC13D)

**Species/Host** Mouse

**Reactivity** Human, Mouse, Rat

**Conjugation** Unconjugated

**Tested Applications** WB

**Immunogen** Full length human recombinant protein of human UNC13D(NP\_954712) produced in HEK293T cell.

**Target** UNC13D

**Note** For research use only

**Application notes** For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

**Isotype** IgG1

**Clonality** Monoclonal


**MW** 123.1 kDa

**Uniprot ID** [Q70J99](#)

**NCBI** [NM\\_199242](#)

**Dilution Range** WB 1:200~1000

**Expiration Date** 12 months from date of receipt.



Western blot analysis of HEK293T cells w...