

Product Datasheet

NKp46 antibody (orb13333)

Description

NKp46 antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

ICC, IF, IHC-P, WB

Immunogen

KLH conjugated synthetic peptide derived from human NKp46. Please contact us for the exact immunogen sequence. The peptide is available as orb374761.

Target

NKp46

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

34 kDa

Uniprot ID

[Q9Z0H5](#), [O76036](#), [Q8C567](#)

NCBI

[001145457](#), [21](#), [001138929](#)

Entrez


[9437](#)

Dilution Range


WB: 1:100-500, IHC-P: 1:50-500, IF/ICC: 1:50-500

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of mouse lymph node



Western blot analysis of rat lymph node tissue



IHC-P image of rat thymus tissue using NKp46 antibody