


Product Datasheet

HORMAD2 antibody (orb1333316)

Description	Mouse monoclonal antibody to HORMAD2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human HORMAD2(NP_689723) produced in HEK293T cell.
Target	HORMAD2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	35.1 kDa
Uniprot ID	Q8N7B1
NCBI	NM_152510
Dilution Range	WB 1:4000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...