

Product Datasheet

Fibrinogen gamma chain (FGG) antibody (orb1333349)

Description Mouse monoclonal antibody to Fibrinogen gamma chain (FGG)

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications IF, WB

Immunogen Full length human recombinant protein of human FGG(NP_000500) produced in HEK293T cell.

Target FGG

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG1

Clonality Monoclonal


MW 46.4 kDa

Uniprot ID [P02679](#)

NCBI [NM_000509](#)

Dilution Range WB 1:2000, IF 1:100

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...