


Product Datasheet

PD-L1 (CD274) antibody (orb1333381)

Description	Mouse monoclonal antibody to PD-L1 (CD274)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB, xMAP
Immunogen	Full length human recombinant protein of human CD274(NP_054862) produced in HEK293T cell.
Target	CD274
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	31 kDa
Uniprot ID	Q9NZQ7
NCBI	NM_014143
Dilution Range	WB 1:4000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...