


Product Datasheet

MT (MCAT) antibody (orb1333389)

Description	Mouse monoclonal antibody to MT (MCAT)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human MCAT(NP_775738) produced in HEK293T cell
Target	MCAT
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	40.6 kDa
Uniprot ID	Q8IVS2
NCBI	NM_173467
Dilution Range	WB 1:400~4000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...