

Product Datasheet

COX4NB (EMC8) antibody (orb1333393)

Description	Mouse monoclonal antibody to COX4NB (EMC8)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human COX4NB(NP_006058) produced in HEK293T cell.
Target	EMC8
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	23.6 kDa
Uniprot ID	O43402
NCBI	NM_006067
Dilution Range	WB 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...