


Product Datasheet

SERPINB3 antibody (orb1333433)

Description	Mouse monoclonal antibody to SERPINB3
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human SERPINB3(NP_008850) produced in HEK293T cell.
Target	SERPINB3
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	44.4 kDa
Uniprot ID	P29508
NCBI	NM_006919
Dilution Range	WB 1:4000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...