



Product Datasheet

NKp30 antibody (orb13335)

Description	NKp30 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human CD337 (81-180/201 aa)
Target	NCR3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	20 kDa
Uniprot ID	O14931
Dilution Range	Flow-Cyt=2ug/Test
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



IHC-P of human cancer tissue (1:600 dilu...