


Product Datasheet

UPRT antibody (orb1333508)

Description	Mouse monoclonal antibody to UPRT
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human UPRT(NP_659489) produced in HEK293T cell.
Target	UPRT
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	33.6 kDa
Uniprot ID	Q96BW1
NCBI	NM_145052
Dilution Range	WB 1:4000, IHC 1:150, IF: 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...