


Product Datasheet

p60 CAF1 (CHAF1B) antibody (orb1333531)

Description	Mouse monoclonal antibody to p60 CAF1 (CHAF1B)
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human CHAF1B(NP_005432) produced in HEK293T cell.
Target	CHAF1B
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	61.3 kDa
Uniprot ID	Q13112
NCBI	NM_005441
Dilution Range	WB 1:200~1000, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.



immunofluorescent staining of COS7 cells...