


Product Datasheet

SIVA (SIVA1) antibody (orb1333542)

Description	Mouse monoclonal antibody to SIVA (SIVA1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human SIVA1(NP_006418) produced in HEK293T cell.
Target	SIVA1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	18.5 kDa
Uniprot ID	O15304
NCBI	NM_006427
Dilution Range	WB 1:1000, IF 1:100, IHC: 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...