

Product Datasheet

MMP13 antibody (orb1333554)

Description	Mouse monoclonal antibody to MMP13
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 104-471 of human MMP13(NP_002418) produced in HEK293T cell.
Target	MMP13
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	42.2 kDa
Uniprot ID	P45452
NCBI	NM_002427
Dilution Range	WB 1:2000, IF 1:100
Expiration Date	12 months from date of receipt.

