


Product Datasheet

Placental lactogen (CSH1) antibody (orb1333562)

Description	Mouse monoclonal antibody to Placental lactogen (CSH1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human CSH1(NP_072166) produced in HEK293T cell.
Target	CSH1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	28.6 kDa
Uniprot ID	P0DML2
NCBI	NM_022640
Dilution Range	WB 1:4000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...