


Product Datasheet

TUBA6 (TUBA1C) antibody (orb1333608)

Description	Mouse monoclonal antibody to TUBA6 (TUBA1C)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human TUBA1C(NP_116093) produced in HEK293T cell.
Target	TUBA1C
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	49.7 kDa
Uniprot ID	Q9BQE3
NCBI	NM_032704
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...