


Product Datasheet

EPCAM antibody (orb1333628)

Description	Mouse monoclonal antibody to EPCAM
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human EPCAM(NP_002345) produced in HEK293T cell.
Target	EPCAM
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	32.7 kDa
Uniprot ID	P16422
NCBI	NM_002354
Dilution Range	WB: 1:1000-1:4000, IF 1:100, Flow 1:100, IHC: 1:150-1:500
Expiration Date	12 months from date of receipt.



immunofluorescent staining of HeLa cells...