


Product Datasheet

Tyrosine Hydroxylase (TH) antibody (orb1333660)

Description	Mouse monoclonal antibody to Tyrosine Hydroxylase (TH)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human TH(NP_000351) produced in HEK293T cell.
Target	TH
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	P07101
NCBI	NM_000360
Dilution Range	WB 1:4000, IF 1:100, IHC: 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...