


Product Datasheet

IL1 beta (IL1B) antibody (orb1333690)

Description	Mouse monoclonal antibody to IL1 beta (IL1B)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human IL1B(NP_000567) produced in HEK293T cell.
Target	IL1B
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	30.6 kDa
Uniprot ID	P01584
NCBI	NM_000576
Dilution Range	WB 1:1000~4000, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...