


Product Datasheet

CD64 (FCGR1A) antibody (orb1333739)

Description	Mouse monoclonal antibody to CD64 (FCGR1A)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human FCGR1A(NP_000557) produced in HEK293T cell.
Target	FCGR1A
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	42.5 kDa
Uniprot ID	P12314
NCBI	NM_000566
Dilution Range	WB 1:4000, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...