


Product Datasheet

CD33 antibody (orb1333746)

Description	Mouse monoclonal antibody to CD33
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human CD33(NP_001763) produced in HEK293T cell.
Target	CD33
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	38 kDa
Uniprot ID	P20138
NCBI	NM_001772
Dilution Range	WB 1:4000, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...