


Product Datasheet

Cytochrome C Oxidase subunit VIc (COX6C) antibody (orb1333790)

Description	Mouse monoclonal antibody to Cytochrome C Oxidase subunit VIc (COX6C)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human COX6C(NP_004365) produced in HEK293T cell.
Target	COX6C
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	8.6 kDa
Uniprot ID	P09669
NCBI	NM_004374
Dilution Range	WB 1:4000, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...