

Product Datasheet

Chromogranin A (CHGA) antibody (orb1333856)

Description	Mouse monoclonal antibody to Chromogranin A (CHGA)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human CHGA(NP_001266) produced in HEK293T cell.
Target	CHGA
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	48.9 kDa
Uniprot ID	P10645
NCBI	NM_001275
Dilution Range	WB 1:2000, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.

