

Product Datasheet

CD45 (PTPRC) antibody (orb1333894)

Description	Mouse monoclonal antibody to CD45 (PTPRC)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human PTPRC(NP_002829) produced in HEK293T cell.
Target	PTPRC
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	147.1 kDa
Uniprot ID	P08575
NCBI	NM_002838
Dilution Range	WB 1:400~4000, IHC 1:500, IF 1:100, FLOW 1:50
Expiration Date	12 months from date of receipt.

