


Product Datasheet

Junctional Adhesion Molecule 1 (F11R) antibody (orb1333901)

Description	Mouse monoclonal antibody to Junctional Adhesion Molecule 1 (F11R)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human F11R(NP_653087) produced in HEK293T cell.
Target	F11R
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	29.8 kDa
Uniprot ID	Q9Y624
NCBI	NM_144504
Dilution Range	WB 1:4000, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...