


Product Datasheet

NY-ESO-1 (CTAG1B) antibody (orb1333957)

Description	Mouse monoclonal antibody to NY-ESO-1 (CTAG1B)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human CTAG1B(NP_001318) produced in HEK293T cell.
Target	CTAG1B
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	17.8 kDa
Uniprot ID	P78358
NCBI	NM_001327
Dilution Range	WB 1:200 - 1:1000, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...