

Product Datasheet

NY-ESO-1 (CTAG1B) antibody (orb1333970)

Description Mouse monoclonal antibody to NY-ESO-1 (CTAG1B)

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications FC, IF, WB

Immunogen Full length human recombinant protein of human CTAG1B(NP_001318) produced in HEK293T cell.

Target CTAG1B

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG1

Clonality Monoclonal


MW 17.8 kDa

Uniprot ID [P78358](#)

NCBI [NM_001327](#)

Dilution Range WB 1:200 - 1:1000, IF 1:100

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...