


Product Datasheet

TRIM45 antibody (orb1333999)

Description	Mouse monoclonal antibody to TRIM45
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 288-580 of human TRIM45(NP_079464) produced in E.coli.
Target	TRIM45
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	64.2 kDa
Uniprot ID	Q9H8W5
NCBI	NM_025188
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of E.coli lysate (...)