


Product Datasheet

RAB21 antibody (orb1334097)

Description	Mouse monoclonal antibody to RAB21
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human RAB21(NP_055814) produced in HEK293T cell.
Target	RAB21
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	24.2 kDa
Uniprot ID	Q9UL25
NCBI	NM_014999
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...