


Product Datasheet

CD272 (BTLA) antibody (orb1334170)

Description	Mouse monoclonal antibody to CD272 (BTLA)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human BTLA(NP_861445) produced in HEK293T cell.
Target	BTLA
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	32.6 kDa
Uniprot ID	Q7Z6A9
NCBI	NM_181780
Dilution Range	WB 1:4000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...