


Product Datasheet

RAB24 antibody (orb1334204)

Description	Mouse monoclonal antibody to RAB24
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human RAB24 (NP_001026847) produced in HEK293T cell.
Target	RAB24
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	22.9 kDa
Uniprot ID	Q969Q5
NCBI	NM_001031677
Dilution Range	WB 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...