


Product Datasheet

ZFYVE21 antibody (orb1334237)

Description	Mouse monoclonal antibody to ZFYVE21
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human ZFYVE21 (NP_076976) produced in HEK293T cell.
Target	ZFYVE21
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	26.3 kDa
Uniprot ID	Q9BQ24
NCBI	NM_024071
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...