

Product Datasheet

C16orf9 (FAM234A) antibody (orb1334267)

Description Mouse monoclonal antibody to C16orf9 (FAM234A)

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications IHC, WB

Immunogen Full length human recombinant protein of human ITFG3 (NP_114428) produced in HEK293T cell.

Target FAM234A

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG1

Clonality Monoclonal


MW 59.5 kDa

Uniprot ID [Q9H0X4](#)

NCBI [NM_032039](#)

Dilution Range WB 1:2000, IHC 1:150

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...