


Product Datasheet

CD272 (BTLA) antibody (orb1334287)

Description	Mouse monoclonal antibody to CD272 (BTLA)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human BTLA(NP_861445) produced in HEK293T cell.
Target	BTLA
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	32.6 kDa
Uniprot ID	Q7Z6A9
NCBI	NM_181780
Dilution Range	WB 1:2000, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...