


Product Datasheet

PSMA2 antibody (orb1334295)

Description	Mouse monoclonal antibody to PSMA2
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 50-234 of human PSMA2(NP_002778) produced in E.coli.
Target	PSMA2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	25.7 kDa
Uniprot ID	P25787
NCBI	NM_002787
Dilution Range	WB 1:2000, IHC 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts (35ug)...