


Product Datasheet

GADD45G antibody (orb1334337)

Description	Mouse monoclonal antibody to GADD45G
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human GADD45G(NP_006696) produced in HEK293T cell.
Target	GADD45G
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	16.9 kDa
Uniprot ID	O95257
NCBI	NM_006705
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...