


Product Datasheet

PCDH7 antibody (orb1334381)

Description	Mouse monoclonal antibody to PCDH7
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human PCDH7(NP_115833) produced in HEK293T cell.
Target	PCDH7
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	132.4 kDa
Uniprot ID	O60245
NCBI	NM_032457
Dilution Range	WB 1:500, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...