


Product Datasheet

CD1C antibody (orb1334408)

Description	Mouse monoclonal antibody to CD1C
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human CD1C(NP_001756) produced in HEK293T cell.
Target	CD1C
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	35.7 kDa
Uniprot ID	P29017
NCBI	NM_001765
Dilution Range	WB 1:500~2000, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...