

Product Datasheet

PPIL3 antibody (orb1334413)

Description	Mouse monoclonal antibody to PPIL3
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human PPIL3(NP_572028) produced in E.coli.
Target	PPIL3
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	18 kDa
Uniprot ID	Q9H2H8
NCBI	NM_131916
Dilution Range	WB 1:2000, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...