


Product Datasheet

CESK1 (CCT8L2) antibody (orb1334471)

Description	Mouse monoclonal antibody to CESK1 (CCT8L2)
Species/Host	Mouse
Reactivity	Human, Monkey
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human CCT8L2(NP_055221) produced in HEK293T cell.
Target	CCT8L2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	59.2 kDa
Uniprot ID	Q96SF2
NCBI	NM_014406
Dilution Range	WB 1:200~2000, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...