


Product Datasheet

Mitochondrial ribosomal protein L11 (MRPL11) antibody (orb1334474)

Description	Mouse monoclonal antibody to Mitochondrial ribosomal protein L11 (MRPL11)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human MRPL11(NP_057134) produced in E.coli.
Target	MRPL11
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	20.5 kDa
Uniprot ID	Q9Y3B7
NCBI	NM_016050
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...