


Product Datasheet

MRPS34 antibody (orb1334516)

Description	Mouse monoclonal antibody to MRPS34
Species/Host	Mouse
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human MRPS34(NP_076425) produced in E.coli.
Target	MRPS34
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	25.5 kDa
Uniprot ID	P82930
NCBI	NM_023936
Dilution Range	WB 1:500~1000, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...