


Product Datasheet

PI 3 Kinase catalytic subunit gamma (PIK3CG) antibody (orb1334522)

Description	Mouse monoclonal antibody to PI 3 Kinase catalytic subunit gamma (PIK3CG)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human PIK3CG(NP_002640) produced in HEK293T cell.
Target	PIK3CG
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	126.3 kDa
Uniprot ID	P48736
NCBI	NM_002649
Dilution Range	WB 1:2000, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...