


Product Datasheet

p15 INK4b (CDKN2B) antibody (orb1334538)

Description	Mouse monoclonal antibody to p15 INK4b (CDKN2B)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human CDKN2B(NP_004927) produced in HEK293T cell.
Target	CDKN2B
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	14.5 kDa
Uniprot ID	P42772
NCBI	NM_004936
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...