


Product Datasheet

POLR1H antibody (orb1334624)

Description	Mouse monoclonal antibody to POLR1H
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human ZNRD1(NP_740753) produced in HEK293T cell.
Target	POLR1H
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	13.7 kDa
Uniprot ID	Q9P1U0
NCBI	NM_170783
Dilution Range	WB 1:1000, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...