


Product Datasheet

LEPRE1 (P3H1) antibody (orb1334627)

Description	Mouse monoclonal antibody to LEPRE1 (P3H1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human LEPRE1(NP_071751) produced in HEK293T cell.
Target	P3H1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	83.2 kDa
Uniprot ID	Q32P28
NCBI	NM_022356
Dilution Range	WB 1:2000, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...