


Product Datasheet

TSC22D1 antibody (orb1334651)

Description	Mouse monoclonal antibody to TSC22D1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human TSC22D1(NP_006013) produced in HEK293T cell.
Target	TSC22D1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	15.5 kDa
Uniprot ID	Q15714
NCBI	NM_006022
Dilution Range	WB 1:2000, IHC 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...